

SILICON OIL

aadpri

Composition: Silicon oil USP.....10 ml (Polydimethy Siloxae)

Do not autoclave

Do not shake

Contains no preservative

Store at room temperature

Read instruction leaflet before use.

Single Use Only

10 ml

Sterile liquid Syringe

Vial

1000 cSt

For intra ocular use only. Not for I.M. / I.V. Injection

Warning: To be used by registered medical practitioner only.

DESCRIPTION

:

SILICONE OIL is fractionated, purified, sterile and apyrogenic silicone oil designed for prolonged tamponade after surgical treatment for severe retinal detachment.

ADMIXTURE INCOMPATIBILITY: Do not admix with any other substance prior to injection

STORAGE : Store between 2°C and 35°C.

SUPPLY : 10 ml vial

Polydimethylsiloxane 1000 Cs

aadSil

STERILE SOLUTION

STERILE



10 ml Vial

10 ml.

STERILE EO



AADPRI OPHTHALMIC COMPANY

G-15, Mahavir krupa complex, Near C.P. Nagar Cross road,
Bhuyangdev char rasta, Ahmedabad- 380061, Gujarat, India.
Mobile no: +91 9374862303,
Email: aadpri@gmail.com Web: www.aadpri.com

aadpri

Low level of catalyst impurities minimize the risks of interactions with the ocular environment.

Quality product with a proven track record.

Greater stability

Subject to all physical and chemical analysis(H-NMR spectral test, gel permeation chromatography study and gas chromatography) to ensure purity level.