

## SODIUM HYALURONATE OPHTHALMIC SOLUTION

aadPri



**PRE-FILLED SYRINGE 1.0ML.**

Sterile Prefilled Syringe with Sterile single use cannula

# D-HYAL 3%

STERILE SOLUTION



### COMPOSITION:

Sodium Hyaluronate B.P.....3% w/v

Aqueous Buffered vehicle.....q.s.

### Description

D-HYAL 3% (Sodium Hyaluronate 3%) is Preservative free Cohesive Ophthalmic Viscosurgical Device.

### Product Features:

D-HYAL 3% offers Excellent Space Maintenance, Crystal Clear Visibility, Ensure Controlled Capsulorhexis, Superior Corneal Endothelial, Tissue Protection and Optimum Viscosity to offer Resistance to shear stress. High Molecular Weight and Long Chain ensure smooth Surgical Procedure and Visco Clearance

### Indications:

Cataract Surgery, Corneal Transplant Surgery, Glaucoma Filtering Surgery, Retinal Surgery & Vitrectomy

### Dosage:

As directed by Physician

### Contraindications:

No known contraindications to the use of D-HYAL 3% when used as recommended.

### Supply:

Available in Sterile 1 ml Pre-Filled Glass Syringes with 27 G cannula.

### Storage:

For Better performance store between 2° - 30° Degree C. D-HYAL 3% should be kept at room temperature for 15 minutes before use.



aadPri

### AADPRI OPHTHALMIC COMPANY

G-15, Mahavir krupa complex, Near C.P. Nagar Cross road,  
Bhuyangdev char rasta, Ahmedabad- 380061, Gujarat, India.

Mobile no: +91 9374862303,

Email: aadpri@gmail.com Web: [www.aadpri.com](http://www.aadpri.com)