


R_x

HYALURONIDASE INJECTION I.P. (OVINE)



COMPOSITION:

10 X 1500 I.U. VIALS

Each vial contains
Hyaluronidase I.P.(Ovine) 1500 I.U.
Prior to use dissolve the lyophilized
power in 1ml of Sterile water for injection

Use immediately after reconstitution,
Do not store.

DOSAGE: As directed by the physician

Not for I.V. Use



AADPRI OPHTHALMIC COMPANY

G-15, Mahavir krupa complex, Near C.P. Nagar Cross road,
Bhuyangdev char rasta, Ahmedabad- 380061, Gujarat, India.
Mobile no: +91 9374862303,
Email: aadpri@gmail.com Web: www.aadpri.com

aadPRI

